



**User manual**  
**for**  
**Child-Pugh Score**

Version **2**, **December 2023**, in English

## 1. The Evidencio platform

The Evidencio platform facilitates the creation, use, validation and implementation of medical prediction models and clinical decision support tools. This user manual specifically relates to the **Child-Pugh Score**. The User Manual can also be referred to as the Instructions For Use (IFU).

Throughout this manual CE-marked content and the term medical device are used interchangeably.

## 2. Disclaimer

Evidencio provides information, models, calculators, equations, and algorithms (tools) intended for use by healthcare professionals. Some of these tools have been certified as CE-medical devices. For such CE-marked content the 'Official Legal Disclaimer for CE-marked content' applies. All other content and tools provided by Evidencio are explicitly only covered by the 'Official Legal Disclaimer for non CE-marked content' both are available here: <https://www.evidencio.com/disclaimer>

## 3. Warnings



### 1. Warnings for CE-marked content

Calculations alone should never dictate patient care, and are no substitute for professional judgement. This tool is only to be used by **physicians** in a clinical setting, and is not for patient use.

Always read the intended use before using this tool.


Before reading the result, double check the filled in values to prevent errors.

Results that concern risk percentages, do not guarantee certain outcomes. When there is a risk present, do not expect an event to not occur at all, even if the risk is very small.

This model is only intended for use in settings where the usage and result of a model are never immediately needed.

## 4. Model landing page

The medical device model on the Evidencio platform is shown in Figure 1. The model landing page contains the following sections, that are indicated in Figure 1.



**A.** Combined Model for End-stage Liver Disease score (Combined MELD Score): (Original MELD, MELD Na (UNOS/OPTN) and MELD 3.0)

**B.** The Model for End-Stage Liver Disease (MELD) scores are a set of scoring tools that quantify end-stage risk disease for transplant planning. The Original MELD and MELD Na (UNOS/OPTN) predict the risk of 3-month mortality. The MELD Na (UNOS/OPTN) includes sodium as input as opposed to the original MELD score. Moreover, this iteration was developed and is supported by the United Network for Organ Sharing (UNOS) and the Organ Procurement and Transplantation Network (OPTN). The MELD 3.0 score predicts the 90-day mortality risk of patients.

**C.** Research authors: W Ray Kim, Ajitha Mannalithara, Julie K Heimbach, Patrick S Kamath, Sumeet K Asrani, Scott W Biggins, Nicholas L Wood, Sommer E Gentry, Allison J Kwong  
Version: 1.0

**D.** Draft | Hepatology | R-Script model

**E.** LOT V-1.0-10085.23.12.20 CE

**F.** UDI (01)008720938015038(8012)v1.0(4326)231220(240)10085

**K.** Download the User manual for Medical device prediction models and consult the Intended use.

**L.** MELD version: MELD 1.0 | MELD Na (UNOS/OPTN) | MELD 3.0

**M.** ... points.

Set all parameters to calculate prediction.

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer.

**G.** Details  
**H.** Intended use  
**I.** Electronic label  
**J.** Release Notes  
**K.** User manual

Figure 1. An example of a model landing page.

### A. Model title

This is the title and name of the model.

### B. Model description

This is a short description of the model.

### C. Research authors

These are the research authors of the paper that originally published the model.

## D. Model tags

These are the tags that are assigned to the model. Evidencio has the following status tags: “Draft”, “Public”, “Private”, “Under review”. Evidencio has the following model type tags: “Composite model”, “Sequential model”, “API model”. Evidencio has the following calculation method tags: “Linear model”, “Logistic regression”, “Cox regression”, “RScript” and “Custom model”. Next to this, there are tags that indicate the specialty e.g. “Cardiology”.

## E. LOT number

The LOT number indicated the model version, the model identifier, and the model publication date. Publication date is indicated as YY.MM.DD.

Additionally, the CE mark is displayed next to the LOT number. This way, medical devices can be easily recognized.

## F. UDI number

The UDI number is an international tool that helps users identify and find information on products. UDI stands for Unique Device Identifier. Evidencio’s UDIs have the following format:

*(01)UDI-DI number(8012)versionnumber(4326)releasedate(240)identificationnumber*

The UDI-DI number is a unique numeric code. For each medical device of Evidencio, a unique UDI-DI is ascribed. This UDI-DI is used as an “access key” for information stored in a unique device identification database (UDID). Information on Evidencio’s medical devices can be found by searching for the UDI-DI number in the following data base: <https://gepir.gs1.org/index.php/search-by-gtin>

## G. Details button

On the top right of the model page, several clickable buttons are displayed that show a pop-up when clicked. The first button opens a pop-up concerning additional information about the model. This pop-up has three sections: Details, Study characteristics and Supporting publications & related files.

### Details

The first part of the additional information concerns the details of the model as shown in Figure 2.



Details		Status	Draft	
Model author	T. Pasma	Status	Draft	
Model ID	10085	Share		
Version	1.0			
Revision date	2023-12-20			
Specialty	Hepatology			
Model type	R-Script model (Calculation)			
MeSH terms	<ul style="list-style-type: none"> <li>• End Stage Liver Disease</li> <li>• Hepatology</li> <li>• Creatinine</li> <li>• Sodium</li> <li>• INR</li> <li>• Bilirubin</li> <li>• Risk</li> <li>• Mortality</li> <li>• Liver</li> <li>• Dialysis</li> </ul>			

Figure 2. The model details.

### Study characteristics

Below the ‘Details section’ the section labeled ‘Study characteristics’ provides information on the characteristics of the patient data used to derive and validate the model. Additional information is provided on the methods used to develop and/or validate the model.

An important part of the Study characteristics is the information on Supporting publications and related files. These sections can be found at the bottom of the Details-pop-up as shown in Figure 3.

## Supporting Publications

<p><b>Title or description</b>          (Development paper MELD 3.0 Score) MELD 3.0: The Model for End-Stage Liver Disease Updated for the Modern Era          DOI: 10.1053/j.gastro.2021.08.050          The new liver allocation system: moving toward evidence-based transplantation policy          DOI: 10.1053/jlts.2002.35927</p>	<p><b>Tags</b>          Model updating          Model updating</p>
--	--

## Related files




Preview	Name	Tags
	Original MELD - A Model to Predict Survival in Patients With End-Stage Liver Disease.pdf 160.59 kB	Original calculator
	MELD UNOS&OPTN - The New Liver Allocation System Moving Toward Evidence-Based Transplantation Policy.pdf 1.43 MB	Model updating
	MELD 3.0 - The Model for End-stage Liver Disease Updated for the Modern Era.pdf 480.35 kB	Model updating

Figure 3. An example of Supporting publications & related files.

## H. Intended use button

### Intended medical use

The **Child-Pugh Score** is intended to be used by professional users who are capable of operating the device and interpreting its results. It can be used to estimate the probability of mortality in patients with end-stage liver disease to support prediction of the 3-month mortality or 90-day survival prognosis of patients.

The three different MELD scores within the Child-Pugh Score largely overlap in their input and outcomes, but do exhibit differences. The MELD 1.0 score combines the values of bilirubin, creatinine and international normalized ratio (INR) to predict the 3-month mortality risk of a patient. The MELD Na (UNOS/OPTN) largely has the same inputs and the outcome is the same, but it uses the result of the MELD 1.0 score in a formula together with sodium-values as an additional variable. The MELD 3.0 score further adds sex and albumin values as variables, and calculates the 90-day survival chance instead of the 3-month mortality risk.

The device is intended to be used for patients with end-stage liver disease by professionals in a clinical setting. The device is not intended for use by patients on their own.

The **Child-Pugh Score** is not intended to replace clinical decision-making, it can only provide information to the user on the mortality risk or survival chance. The user can use this information to support clinical decision-making regarding prognosis and treatment of the patient. In practice, treatment of the patient usually entails palliative care or liver transplantation.

### Clinical Benefit

The **Child-Pugh Score** is intended to assist patients with relevant and specified clinical outcome parameters. Concretely, this is achieved by estimating a risk in order to support clinical decision-making aimed at end-stage liver disease patients, in order to support clinical decision-making regarding patient prognosis. Correct functioning of the **Child-Pugh Score** can result in these clinical benefits:

- The **Child-Pugh Score** can aid a professional in providing a patient with an accurate prognosis. More accurate prognosis can support necessary decision-making of patients with end-stage liver disease and organization of their life, such as the need for palliative care.
- Use of the **Child-Pugh Score** can positively impact patient management by supporting the decision-making concerning application for a liver transplantation.
- [Adapt, add or remove clinical benefits depending on those of the SaMD in question]
- Digital implementation of the algorithm underlying the **Child-Pugh Score** as a medical device can improve the speed and reliability of calculation. This would further increase the accuracy of the prognosis and by extent increase the chance for the above-mentioned benefits.

## Intended target population and exclusion

The **Child-Pugh Score** should be used for **end-stage liver disease patients**. The **Child-Pugh Score** is intended to be used only for a specific group of patients, corresponding to the below indications and contra-indications. The target population of the model is **end-stage liver disease patients**, provided that they fit the listed indications and contra-indications.

## Clinical indication

The **Child-Pugh Score** should be used for patients who meet the following inclusion criteria:

- Patients should be at least 18 years or older \*
- Having end-stage liver disease

*\*The authors of the MELD score suggest that it can be used for patients of 12 years or older. However, this is not universally accepted, e.g. the OPTN policy acknowledges three groups: less than 12 years, 12 or older and 18 or older. In contrast, Eurotransplant uses an age limit of 18 years or older, below which the PELD is recommended. Because of this lack of consensus, the age for the use of the Child-Pugh Score is set to 18 years or older.*

## Contra-indications

The **Child-Pugh Score** should not be used for patients who meet one or more of the following exclusion criteria:

- Patients younger than 18 years (See also the inclusion criteria)

Besides, there is a large discussion regarding certain medical conditions, known as 'MELD Exceptions'. Many liver allocation and transplantation programs provide additional points to patients with these conditions to compensate for them for the purpose of equity among patients regarding liver allocation. This measure is due to the fact that the MELD score results are considered to not accurately represent the mortality risk of patients with these conditions. However, there is no consensus on the exact list of conditions that should be regarded as exceptions as well as how to deal with them (e.g. the amount of points to be distributed) to reach equity for all patients.

The following conditions are acknowledged as MELD Exceptions in different liver allocation and transplantation programs, and thus the use of the Child-Pugh Score for patients with these conditions should be performed with caution and the MELD Score results should be interpreted within the context of the condition:

Generally acknowledged MELD Exceptions:

- Hepatopulmonary syndrome (HPS)
- Portopulmonary hypertension (PPH)
- Primary hyperoxaluria
- Cystic fibrosis (CF)
- Cholangiocarcinoma (CCA)
- Hepatocellular carcinoma (HCC) (often together with Milan criteria to determine if HCC patients are eligible/suitable for transplantation)

Additionally, the ELITA acknowledges the following exceptions:

- Cholangitis (Primary sclerosing cholangitis and biliary sepsis/secondary sclerosing cholangitis)
- Neuroendocrine tumors
- Polycystic liver disease (PLD)
- Hepatic artery thrombosis (HAT)
- Persistent hepatic dysfunction (incl. small-for-size syndrome)
- Hereditary hemorrhagic teleangiectasia (Rendu-Osler-Weber-Syndrome)
- Biliary atresia
- Non-metastatic hepatoblastoma
- Urea-cycle disorder/organic acidemia
- Hepatic hemangioendothelioma

Moreover, the following exceptions are mentioned in other literature:

- Familial amyloid polyneuropathy (FAP)
- (Unusual) metabolic disease
- Hepatorenal syndrome (HRS) (Eurotransplant region)
- Amyloidosis

Acknowledged by the 2006 MELD Exception Study Group and Conference (MESSAGE):

- Unusual tumor
- Unusual metabolic disease
- Ascites

- Hepatic encephalopathy
- Gastrointestinal bleeding
- Budd-Chiari Syndrome (BCS)
- Pruritis

Moreover, in an OPTN proposal for addition of serum sodium, it is stated that the serum sodium values for patients with hyperglycemia should be corrected via an additional formula. Thus, caution is also necessary if the Child-Pugh Score is used for these patients.

### User profile

Since End-stage liver disease and mortality are both regarded as a 'critical healthcare situation or condition', the use of the SaMD is intended for specialised trained users in a clinical setting. The Child-Pugh Score should not be used by patients alone. Users do not require additional training prior to the use of the medical device.

### Intended use environment

The SaMD can be used as made available on the Evidencio platform in any actively supported web-browser on personal computers, mobile devices, or tablet PCs, and on the mobile app provided by Evidencio. Furthermore, the SaMD can be used through the Evidencio iFrame representation of the SaMD, as an embedded view, provided that the specific Evidencio guidelines for iFrame implementations of this SaMD are adhered to. In addition, automated calculation of the combined MELD score is allowed through Evidencio's API. The model is only intended for use in settings where the usage and result of a model are never immediately needed.

### Functioning, physical principle

The underlying models that comprise the Child-Pugh Score concern linear models, while Cox proportional-hazards regression models are used for the conversion from the score result to a risk percentage. The MELD score algorithms and the accompanying conversion algorithms are combined in a single R-script. The acquisition and processing of the data, the analyses to assemble the relevant criteria for the MELD 1.0, MELD Na (UNOS/OPTN) and MELD 3.0 scores, as well as their setup and refinement, are described in the original studies from Malinchoc, Kamath, Freeman and Wiesner et al. (MELD 1.0), Biggins et al., Kim et al., the 2013 proposal of the Liver and Intestinal Organ Transplantation Committee, Alcorn et al., Sharma et al., and Kalra et al. (MELD Na (UNOS/OPTN)), and Kim et al. (MELD 3.0), respectively.

Entering the details of an individual in the web-application of Evidencio initiates the calculation of the model and provides the total number of points and the corresponding risk percentage of the patient.

## I. Electronic label button

The electronic label button opens a pop-up with the location and address of Evidencio, the LOT number, the UDI number, the CE-mark, the medical device logo and a download link for the declaration of conformity of the medical device. The example of the electronic label is shown in Figure 4.




Figure 4. Example of the electronic label

## J. Release notes

The 'Release Notes' button opens a pop-up with the latest release notes of the model. Here you can find what has changed over the last versions of the model. Additionally, if there are any known residual anomalies the user should be aware of, they are listed here.

## K. User Manual

This user manual can be found in three places: 1) under the short description, 2) on the right of the model page, and 3) in the electronic label. Additionally, all versions of the user manual can be found in the general page for all user manuals for medical devices. The page can be found under the 'About' drop-down menu button as shown in Figure 5. The user manual page is shown in Figure 6.




Figure 5. The drop-down menu where the user manual page can be found.




Figure 6. The user manual page for all user manuals.

You (The user) can always print this downloaded manual. When necessary, you can request a paper version of the manual to be sent to you by mail. Evidencio's contact details are listed in Chapter 6 of this user manual.

## L. Input section

The Evidencio platform allows two separate input variables; categorical, and continuous variables.

### Categorical variables

In the example shown in Figures 7 and 8, the **MELD version** variable concerns a categorical variable. The version that is wished to be used can be entered by clicking on either button. The selected button changes to green, as seen in Figure 8.




Figure 7. The variable for MELD version, where no button has been clicked, and thus no input has been provided by the user.




Figure 8. The variable for MELD version, where the "MELD 3.0" button has been clicked.



## Continuous variables

In the example shown in Figure 9, the **Serum bilirubin** variable, exemplifies a **continuous variable**. The plausible ranges for the variables are used for the model.

The details for a patient can be entered by sliding the button to the correct value, or by entering the correct value in the box on the right-hand side (i.e., where the **20** is entered for **Serum bilirubin**).



**Serum bilirubin**

Values below 0.1 have to be set as 1.0


1 100

20 mg/dL

Figure 9. The variable for Serum bilirubin, where "20" has been entered

## Details on variable measurements

Directly underneath the name for each variable, additional details can be provided on the methods required to enter the correct value for each variable. In Figure 10, the details below **Serum sodium** explain what the ranges of the variable is.



**Serum sodium**

Sodium is limited to a range of 125-137 mmol/L, and if outside of these bounds, is set to the nearest limit.

125 137

mmol/L

Figure 10. An example on how additional information can be provided for a variable.

## M. Result section


At the bottom of the page, the results of the model are shown.

### Result calculation

When all variables are filled in, a result will be calculated. No risk is displayed until all variables are filled in. The result section indicates "Set all parameters to calculate prediction."

### Result interpretation

In the result interpretation, a risk stratification is given based on the risk score. An example of the information is shown In Figure 11.



**26 points.**

	Result
MELD score	26
MELD (estimated 90-day survival)	85.4 %

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer.

Figure 11. The result information

## Relevant information for correct use of the model

At the bottom of the page, there is a link to Evidencio's terms and conditions of use, the privacy policy, and the disclaimer.

## 5. Use of Medical devices

In general, and unless explicitly stated otherwise, CE-marked tools on Evidencio are only to be used by **physicians** in a clinical setting, and are not for patient use.

To use the tool, Evidencio requires a stable internet connection and runs on the following devices:

- Personal computers or laptops using the following browsers:
  - Safari (the latest three versions)
  - Chrome (the latest three versions)
  - Firefox (the latest three versions)
  - Edge (the latest three versions)
- Tablets or smartphones running on the next operating systems:
  - IOS (the latest three versions)
  - Android (the latest three versions)

The medical device cannot be used in combination with Internet Explorer. The personal computers, laptops, tablets or smartphones used should at least be able to have an internet connection and use the browsers mentioned above. The minimal screen resolution should be 800x600.

Furthermore, the model may be used through the Evidencio iFrame representation of the calculator, as an embedded view, provided that the specific Evidencio guidelines for iFrame implementations of that model are adhered to.

The Evidencio SaMD models can be used with any browser settings that don't distort the regular display of websites, with a 50% to 500% zoom rate, and at a display resolution starting from 800x600. However, factory recommended browser settings, 100% zoom rate and regular display resolution are recommended.

This model is only intended for use in settings where the usage and result of a model are never immediately needed.

## 6. Manufacturer details

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the country in which you, the reader, are established. A competent authority is the institute that governs all issues related to medical devices in a country.

Contact details of your competent authority can be found here: <https://www.ema.europa.eu/en/partners-networks/eu-partners/eu-member-states/national-competent-authorities-human>

Please contact Evidencio when you suspect any malfunction or changes in the performance of a medical device. Do not use the device, until Evidencio replies to your message that it is safe to start using it again.

Contact details of Evidencio:



Evidencio B.V., Irenesingel 19, 7481 GJ Haaksbergen, The Netherlands  
www.evidencio.com  
tel: +31 53 85195 08  
e-mail: info@evidencio.com