



# User manual for the Child-Pugh Score

Version 2.0, JUN 2026, in English



## Table of Contents

- 1. The Evidencio platform ..... 3
- 2. Disclaimer ..... 3
- 3. Warnings for CE-marked content..... 3
  - 3.1. Notice to the user..... 4
- 4. Device description Child-Pugh Score..... 4
  - 4.1. Lifetime, residual risks and side effects ..... 4
- 5. Electronic label ..... 5
  - 5.1. LOT number ..... 5
  - 5.2. UDI number ..... 5
- 6. Intended purpose ..... 5
  - 6.1. Intended use..... 5
  - 6.2. Clinical benefit..... 6
  - 6.3. Indented target population and exclusion..... 6
    - 6.3.1. Clinical indication ..... 6
    - 6.3.2. Contra-indication..... 6
  - 6.4. User profile ..... 6
  - 6.5. Intended use environment ..... 6
  - 6.6. Physical interaction..... 6
  - 6.7. History/ versions of the MDSW ..... 7
- 7. Result interpretation..... 7
- 8. Additional information ..... 7
  - 8.1. Details..... 7
  - 8.2. Input variables..... 8
  - 8.3. Algorithm..... 8
  - 8.4. Study characteristics..... 8
  - 8.5. Supporting publication & Related files ..... 9
  - 8.6. Analytical performance characteristics ..... 9
  - 8.7. Clinical performance characteristics..... 9
  - 8.8. Release notes..... 10
- 9. Using the algorithm on the Evidencio website ..... 11
  - 9.1. General algorithm landing page..... 11
- 10. Implementation of the algorithm through an API ..... 19
- 11. User manual revision history ..... 21
- 12. Manufacturer details..... 21

## 1. The Evidencio platform

The Evidencio platform facilitates the creation, use, validation and implementation of medical prediction algorithms and clinical decision support tools. This User Manual specifically relates to the Child-Pugh Score. The User Manual can also be referred to as the Instructions For Use (IFU). The Child-Pugh Score meets the requirements laid down in: Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU. Compliance with the applicable regulations is provided by means of declaration of conformity.

Throughout this manual CE-marked content and the term medical device are used interchangeably.

## 2. Disclaimer

Evidencio provides certain CE-marked information, calculators, equations, and algorithms (tools) on any of its websites, applications, apps, or services. These tools may only be used in accordance with the intended use / intended purpose that has been published with the respective CE-marked tool.

In general, and unless explicitly stated otherwise, CE-marked tools on Evidencio are only to be used by healthcare professionals and are not for patient use.

The CE-marked content on the platform is to be regarded as a specific set of tools, apart from the general platform content. Any available content, on any of the websites, applications, apps, or services provided by Evidencio that is not clearly labelled as a CE-marked tool is explicitly not covered by this disclaimer for CE-marked content, the general Evidencio Disclaimer for non-CE-marked content applies.

CE-marked tools may provide limited professional advice to the intended user(s). However, the intended user must exercise their clinical judgment as to the information these tools provide.

Evidencio does not assume any liability or responsibility for damage or injury (including death) to you, other persons, or property arising from any misuse of any product, information, idea, or instruction contained in the tools provided to you.

The disclaimer for non-CE-marked content is available on the Evidencio website: <https://www.evidencio.com/disclaimer>.

Your use of the websites, applications, apps, or services provided by Evidencio is subject to our Terms & Conditions, which can be found here: <https://www.evidencio.com/terms-conditions>.



## 3. Warnings for CE-marked content

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer on: <https://www.evidencio.com/disclaimer>. This tool is only to be used by healthcare professionals, and is not for patient use.

Always read the intended use before using this tool.

Always make sure the patient complies with the clinical indications and clinical contra-indications as stated on the Evidencio website, and in **paragraphs 6.3.1** and **6.3.2** of this user manual respectively.

Before reading the result, double check the filled in values to prevent errors.

Results that concern risk percentages, do not guarantee certain outcomes. When there is a risk present, do not expect an event to not occur at all, even if the risk is very small. Conversely, a high risk does not guarantee that an event will occur.

This algorithm is only intended for use in settings where the usage and result of an algorithm are never immediately needed.

The data used to perform the calculations is stored by Evidencio to enhance algorithm function and allow issues to be traceable for further improvements. For details, see the privacy policy on our website at: <https://www.evidencio.com/privacy-policy>.

### 3.1. Notice to the user

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the country in which you, the reader, are established. A competent authority is the institute that governs all issues related to medical devices in a country.

Please contact Evidencio when you suspect any malfunction or changes in the performance of a medical device. Do not use the device, until Evidencio replies to your message that it is safe to start using it again.

## 4. Device description Child-Pugh Score

The Child-Pugh Score is an algorithm intended to use Bilirubin, Albumin, Prothrombin time or INR, Encephalopathy grade and Ascites severity to estimate severity of liver disease in patients with liver dysfunction, in order to support healthcare professionals with decisions surrounding liver disease. It is intended for patients with liver disease and not for use in patients younger than 18 years.

The Child-Pugh score was originally developed in 1973. It was then developed to predict mortality during surgery, but is nowadays often used to determine prognosis in patients suffering from chronic liver disease (mostly cirrhosis). The Child-Pugh score is also known as the Child-Turcotte-Pugh score or the Child score.

The MDSW's underlying mathematical formula is a custom calculation. The Child-Pugh Score was developed to predict the severity of liver dysfunction. Use of the algorithm positively impacts patient management by optimising the assessment of the severity of liver dysfunction of patients with liver disease and informing clinical management on further prognostic/therapeutic options.

The calculation of the algorithm is performed by communication with the Evidencio platform, hosted at [www.evidencio.com](http://www.evidencio.com). Local hosting of the algorithm on another platform is possible using Docker images. The algorithm is also accessible by 3rd party applications through the API and iFrame implementation. The Evidencio platform is managed under Evidencio's certified quality management system that ensures the correctness of calculations and availability of its services.

The Summary of Safety and Performance for this device will be made available via EUDAMED once the relevant module is fully operational. In the meantime, the Summary of Safety and Performance can be requested from the manufacturer and will be provided without undue delay.

### 4.1. Lifetime, residual risks and side effects

The Child-Pugh Score is software, and does not expire. The lifetime is initially set at 5 years from certification, if the state of the art does not change in such a way as to negatively affect the benefit-risk of the device, the lifetime can be extended.

No steps are required to be undertaken by the user to decommission a product when it is taken off the market. If the lifetime is not extended, a notice will be placed on the algorithm page on the platform. When a device is taken off the market, users may be informed about this (e.g. through e-mail).

Evidencio has identified a series of risks associated with the use of this algorithm.

The Child-Pugh Score is a low-risk device, there are no noticeable risks involved outside of possible mis-estimation of the severity of liver dysfunction in patients, and all residual risks are accepted.






Most risks can be defined into two main groups, depending on their outcome.

- a) The risk calculation was wrong or;
- b) The MDSW prediction algorithm is inaccessible.

A wrong risk calculation can be the result of erroneous input values or an error in the mathematical calculation. Technical risks, including the erroneous calculations or inaccessibility due to a technical error, have been mitigated when possible. These measures focussed on reducing the risks' probability and severity. Concluding that the risks could not be mitigated further, the residual risks were classified as *low-level and acceptable*. The Child-Pugh Score does not have any direct side effects.

## 5. Electronic label

The electronic label of this device contains the following information:

	<b>Name of the device</b>	Child-Pugh Score
	<b>Manufacturer information</b>	Evidencio B.V., Irenesingel 19, 7481 GJ Haaksbergen, The Netherlands
	<b>LOT number</b>	V-2.2-253.26.06.29
	<b>UDI number</b>	(01)08720938015090(8012)v2.2(4326)260629(240)253
	<b>IVD indication</b>	<i>In vitro</i> diagnostic medical device

The electronic label can be found on the Evidencio website, see also section **I** and **Figura 5** in **Chapter 10**.

The electronic label on the website further contains the option to download the **User Manual** and **Declaration of conformity** (DoC).

### 5.1. LOT number

The LOT number indicates the algorithm version, the algorithm identifier, and the algorithm publication date. Publication date is indicated as YY.MM.DD.

### 5.2. UDI number

Stands for Unique Device Identifier (UDI) number, which is an international tool that helps users identify and find information on products. Evidencio's UDI's have the following format:

*(01)[UDI-DI number](8012)[versionnumber](4326)[releasedate](240)[identificationnumber]*

The UDI-DI (Device Identifier) number is a unique numeric code. For each medical device of Evidencio, a unique UDI-DI is ascribed. This UDI-DI is used as an "access key" for information stored in a unique device identification database (UDID). Information on Evidencio's medical devices can be found by searching for the UDI-DI number in the following data base:

<https://gepir.gs1.org/index.php/search-by-gtin>.

## 6. Intended purpose

### 6.1. Intended use

The Child-Pugh Score is an algorithm intended to use Bilirubin, Albumin, Prothrombin time or INR, Encephalopathy grade and Ascites severity to estimate severity of liver disease in patients with liver dysfunction, in order to support healthcare professionals with decisions surrounding liver disease.

The Child-Pugh Score is medical device software that automates the calculation of the formula. It requires quantitative and qualitative inputs to provide a quantitative output.

The Child-Pugh Score is not intended to replace clinical decision-making; it can only provide information to the healthcare professional on the estimation of the severity of liver dysfunction. The healthcare professional can use this information to

support clinical decision-making regarding prognostic/therapeutic options for the patient. In practice, the Child-Pugh Score is used in a wide variety of specific contexts.

## 6.2. Clinical benefit

The benefits and risks associated with the use of the Child-Pugh score for the patient are indirect. The benefits arise from clinical decisions made using the Child-Pugh score in combination with other clinical and patient-specific factors. The Child-Pugh Score can result in these clinical benefits:

- The Child-Pugh Score can assist in risk stratification based on the severity of liver dysfunction of patients with liver disease and informs clinical management on further prognostic/therapeutic options.
- Use of the algorithm allows pharmaceutical dosage to be adapted to patient liver function as described in the guidance from the U.S. Food and Drug Administration (FDA)\* and the European Medicines Agency (EMA, formerly EMEA)\*\*.

\* U.S. Department of Health and Human Services. (2003). *Pharmacokinetics in Patients with Impaired Hepatic Function: Study Design, Data Analysis, and Impact on Dosing and Labeling*. *Clinical Pharmacology*. \*\* Committee for medicinal products for human use. (2005). *Guideline on the evaluation of the pharmacokinetics of medicinal products in patients with impaired hepatic function*. European Medicines Agency.

## 6.3. Indented target population and exclusion

The Child-Pugh Score is intended to be used for patients who fit the clinical indications and contra-indications listed below.

### 6.3.1. Clinical indication

The Child-Pugh Score should be used for patients who meet the following inclusion criteria:

- Patients should be at least 18 years or older.
- Patients must have a confirmed diagnosis of liver disease leading to cirrhosis.

### 6.3.2. Contra-indication

The Child-Pugh Score should not be used for patients who meet the following exclusion criteria:

- Patients younger than 18 years.

The Child-Pugh Score should not be seen as directly interchangeable with other hepatic function assessment systems such as the NCIc classification.

The Child-Pugh score is designed for chronic liver disease. Its utility in acute liver failure and liver dysfunction primarily due to advanced liver metastases or primary liver malignancies is limited. Other scoring systems may be more appropriate in acute settings.

## 6.4. User profile

The result of the Child-Pugh Score is intended to be reviewed and interpreted by healthcare professionals. Results shall always be reviewed and interpreted by healthcare professionals, in the context of the patient's clinical history and other diagnostic test results. Healthcare professionals do not require additional training prior to the use of the medical device. The device is not intended for use by patients on their own.

## 6.5. Intended use environment

The MDSW can be used as made available on the Evidencio platform in any actively supported web-browser on personal computers, mobile devices, or tablet PCs. Users can manually enter the required input data through the user interface. In addition, the MDSW is available as an embedded view via Evidencio's iFrame representation. Automated calculation of the device is enabled through Evidencio's API. The device is only intended for use in healthcare settings where the immediate application and outcomes of the device are not required. The device is not intended to be used at the bedside of the patient.

## 6.6. Physical interaction

The MDSW is stand-alone software and does not come into contact with any bodily or other material of the patient, user or otherwise.

## 6.7. History/ versions of the MDSW

The version of the Child-Pugh Score concerns the initial version of the MDSW of which Evidencio is the manufacturer.

## 7. Result interpretation

### Primary outcome

The primary output of this device is the Child Pugh Score, which is given as an integer number of points (between 5 and 15). The other important result is the Child-Pugh class (A, B, or C) associated with the score.

### Conditional information

The Child-Pugh Score can be interpreted as the severity of liver dysfunction. **Table 1** shows the classifications corresponding to a certain number of points. The outcome of the Child-Pugh Score can be used to stratify patients according to liver functioning and prognosis. However, the effectiveness of the Child-Pugh Score and the associated risk percentages can differ between different Cirrhosis aetiologies and patient conditions.

**Table 1.** Child-Pugh outcomes.

Points	Classification
5 or 6 points	A
7 to 9 points	B
10 to 15 points	C

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See the Evidencio website for the full disclaimer; <https://www.evidencio.com/disclaimer>.

## 8. Additional information

### 8.1. Details

<b>Algorithm author</b>	Evidencio
<b>Root algorithm ID</b>	253
<b>Version</b>	2.2
<b>Revision date</b>	2026-JUN-29
<b>Speciality</b>	Gastroenterology
<b>Algorithm type</b>	Custom calculation
<b>MeSH terms</b>	<ul style="list-style-type: none"> <li>• Liver Cirrhosis</li> <li>• Chronic Liver Failure</li> <li>• Liver Disease</li> </ul>

## 8.2. Input variables

To perform the calculations successfully, the Child-Pugh Score requires the input variables as listed in **Table 2**. Depending on the type of clotting time measurement chosen either international normalized ratio or prothrombin time is filled in. **Table 3** shows the variables with options/ranges and their corresponding points.

**Table 2.** Variables used as input for the Child-Pugh Score.

Name	Description	Type	Range (step size)	Units
<b>Encephalopathy</b>	Presence of encephalopathy	Categorical	None Grade I-II Grade III-IV	-
<b>Bilirubin</b>	Bilirubin level in $\mu\text{mol/L}$ or $\text{mg/dL}$	Continuous	1.5 – 850 (0.1) 0.1 – 45 (0.1)	$\mu\text{mol/L}$ $\text{mg/dL}$
<b>Albumin</b>	Albumin level in $\text{g/L}$ or $\text{g/dL}$	Continuous	0 – 50 (1) 0 – 5 (0.1)	$\text{g/L}$ $\text{g/dL}$
<b>Clotting time measurement</b>		Categorical	International Normalized Ratio (INR) Prothrombin time	-
<b>INR</b>	International normalized ratio	Continuous	0.5 – 20 (0.1)	ratio
<b>Prothrombin time</b>		Continuous	1 – 100 (0.1)	Seconds
<b>Ascites</b>	Presence of ascites	Categorical	None Mild Moderate to severe	-

**Table 3.** The variables with their corresponding possible values.

Name	Option/range	Points	
<b>Encephalopathy</b>	None	(+1)	
	Grade I – II	(+2)	
	Grade III – IV	(+3)	
<b>Bilirubin</b>	<34.2 $\mu\text{mol/L}$ (<2 $\text{mg/dL}$ )	(+1)	
	34.2 - 51.3 $\mu\text{mol/L}$ (2-3 $\text{mg/dL}$ )	(+2)	
	>51.3 $\mu\text{mol/L}$ (>3 $\text{mg/dL}$ )	(+3)	
<b>Albumin</b>	>35 $\text{g/L}$ (3.5 $\text{g/dL}$ )	(+1)	
	35 - 28 $\text{g/L}$ (3.5 - 2.8 $\text{g/dL}$ )	(+2)	
	<28 $\text{g/L}$ (2.8 $\text{g/dL}$ )	(+3)	
<b>Clotting time measurement</b>	<b>INR</b>	<1.7	(+1)
		1.7 - 2.3	(+2)
		>2.3	(+3)
	<b>Prothrombin time</b>	<4 seconds	(+1)
		$\geq 4 - \leq 6$ seconds	(+2)
		>6 seconds	(+3)
<b>Ascites</b>	None	(+1)	
	Mild	(+2)	
	Moderate to severe	(+3)	

## 8.3. Algorithm

The algorithm of the Child-Pugh Score is as follows:

$$\text{Child Pugh Score} = \text{Encephalopathy} + \text{Bilirubin} + \text{Albumin} + (\text{INR OR Prothrombin time}) + \text{Ascites}$$

**Algorithm of the Child-Pugh Score.**

## 8.4. Study characteristics

The Child-Pugh Score was developed by Pugh *et al.* (1973) with a cohort that consisted of 38 patients seen between January 1966 and July 1972. A systematic review was performed by Cholongitas *et al.* (2005) where the Child-Pugh Score was compared to the end-stage liver disease model, and the Child-Pugh Score is part of drug dosage procedures as described by the FDA (2003) and EMA (2005).

Information on the categorical characteristics of the patient data used to derive the algorithm is provided in **Table 4**.

**Table 4.** Categorical patient characteristics on the derivation cohort.

Name	Subset/Group	Nr. of patients
Child-Pugh Grade	Grade A	7
Child-Pugh Grade	Grade B	13
Child-Pugh Grade	Grade C	18
Hepatic Encephalopathy	Yes	19
Hepatic Encephalopathy	No	19
Jaundice	Yes	19
Jaundice	No	19
Ascites	Yes	17
Ascites	No	21
Mortality	In-hospital death	21
Mortality	Death within 6 months	5
Mortality	Survived for 1 year	9

## 8.5. Supporting publication & Related files

Several relevant studies, such as the original derivation study by Pugh *et al.* (1973) are contained in **Table 5**. These publications have tags to identify their link with the algorithm. Examples of relevant tags are; “Peer review”, “Internal validation”, “External validation”, and “TRIPOD”. Publications that have the tags: “Internal validation” or “External validation”, contain data on the performance characteristics of the device.

**Table 5.** Overview of selection of supporting publications & Related files.

Internal validation Peer review	<p><b>Transection of the oesophagus for bleeding oesophageal varices.</b> (1973) <i>Pugh R.N., Murray-Lyon I.M., Dawson J.L., Pietroni M.C., Williams R.</i></p> <p>DOI: <a href="https://doi.org/10.1002/bjs.1800600817">https://doi.org/10.1002/bjs.1800600817</a></p>
External validation Peer Review	<p><b>Systematic review: the model for end-stage liver disease – should it replace Child-Pugh's classification for assessing prognosis in cirrhosis?</b> (2005) <i>E. Cholongitas, G. V. Papatheodoridis, M. Vangeli, N. Terreni, D. Patch, A. K. Burroughs</i></p> <p>DOI: <a href="https://doi.org/10.1111/j.1365-2036.2005.02691.x">10.1111/j.1365-2036.2005.02691.x</a></p>

## 8.6. Analytical performance characteristics

To demonstrate the analytical performance of the Child-Pugh Score, evidence was collected based on five requirements. This led to the following results:

- A code review and functional test showed that the calculation of the online tool provides the exact same results as described in the paper by Pugh *et al.* (1973).
- Monthly uptime reports show that the device is available online with an uptime of at least 99%.
- The calculation time is within 2 minutes, otherwise an error is given to the manufacturer, this is analysed each 6 months in the analysis of quality data.
- Absence of unacceptable cybersecurity vulnerabilities.
- On a scale from 1 to 5, where 5 means that users trust the implementation very much, the reliability score was 4.44, and the accuracy score was 4.56 on average.

## 8.7. Clinical performance characteristics

The Child-Pugh Score is used for categorisation in a wide variety of patients. It is often used as a proxy for or definition of liver function or decompensation, and to decide upon treatment paths together with other patient characteristics. Therefore, it is difficult to determine performance characteristics. The Child-Pugh Score is often compared with the MELD Score and while they have similar functionality, the Child-Pugh only provides a broad category, whereas the MELD Score is used for a ranking of patients within the list.

### Statistic performance characteristics

In the vast majority of Child-Pugh Score applications, the C-statistic is above 0.6. In the case of those with lower C-statistics, this sometimes coincides with inclusion criteria based in part on Child-Pugh, such as Liver Transplantation. The Child-Pugh Score performs well for patients with among others (Alcohol-related) Cirrhosis, TIPS placement, or Acute-on-Chronic Liver Failure. Kaplan-Meier plots show a clear difference between groups based on Child-Pugh Scores.

### Benefit assessment

The direct impact of the Child-Pugh Score is difficult to quantifiably measure. It is ubiquitous throughout the field of hepatology. The EMA and the FDA include recommendations towards the use of the Child-Pugh Score, especially if no better marker of hepatic functioning is found. There are several alternatives, often performing better when it comes to prognosis, but these are generally more complicated or meant for a more specific sub-group of patients. None of these have displaced the Child-Pugh Score, though the MELD is used in some similar situations. Thus, the main clinical benefit of the Child-Pugh is the ability to compare an individual patient to similar patients within the scientific literature, as well as to provide dosage guidelines for specific pharmacological treatment regimens. Furthermore, the Child-Pugh (and often the different MELD versions) serve a special role as a baseline for comparison with other prognostic tools. As such, their performance is measured in how much better than the Child-Pugh they perform for a group of patients.

Another difficult to quantify benefit of the Child-Pugh Score is that there are other algorithms that have been designed to use the Child-Pugh Score. Additionally, there are models which are designed specifically for use within a specific Child-Pugh Category to add further discriminatory ability.

No cost-effectiveness Analysis has been conducted which includes the Child-Pugh Score.

## 8.8. Release notes

The release notes for each publicly available version of the device can be found on the Evidencio website page for the Child-Pugh Score: <https://www.evidencio.com/models/show/253>, selecting the correct device (version), and clicking on Release Notes. It is recommended to read these notes after a version update to see if these changes are relevant to you. Please make sure the correct algorithm version is selected.

## 9. Using the algorithm on the Evidencio website

Using the tool on the Evidencio website requires a stable internet connection. The tool was developed to work on the four most commonly used internet browsers; Google Chrome (version 135.0.7049.115 and higher), Mozilla Firefox (version 137.0.2 and higher), Microsoft Edge (version 135.0.3179.98 and higher), and Apple Safari (version 18.4 and higher). The medical device cannot be used in combination with Internet Explorer.

The tool can also be accessed on mobile devices running the most recent versions of the Android (version 15 and higher) and iOS (version 18.4.1 and higher) operating systems.

Correct functioning of the tool with earlier versions of these browsers cannot be guaranteed.

The personal computers, laptops, tablets or smartphones used should at least be able to have an internet connection and use the browsers mentioned above.

Furthermore, the algorithm may be used through the Evidencio iFrame representation of the calculator, as an embedded view, provided that the specific Evidencio guidelines for iFrame implementations of that algorithm are adhered to.


The Evidencio MDSW algorithms can be used with any browser settings that don't distort the regular display of websites, with a 50% to 500% zoom rate, and at a display minimal screen resolution starting from 800x600. However, factory recommended browser settings, 100% zoom rate and regular display resolution are recommended.

The MDSW is intended for authorised users only, and should not be used by unauthorised personnel.

This algorithm is only intended for use in settings where the usage and result of an algorithm are never immediately needed.

### 9.1. General algorithm landing page

An example of a medical device algorithm interface on the Evidencio platform is shown in **Figure 1**. The different sections indicated are explained in this chapter.



**A.** Algorithm Name

**B.** Here a description of the algorithm is provided. This will include a short description what the device does, why it was developed, the target population, and the output. It may also mention the current state of the art. Please make sure the patient meets the clinical indications and contra-indications, which can be found in the Intended Purpose section and the User Manual.

**C.** Research authors: First Author Name, Second Author Name, exc. Version: 1.0

**D.** Public | Speciality (e.g. Geriatrics) | Custom calculation | Versions

**E.** LOT V-1.0-10513.25.02.11 CE

**F.** UDI (01)012345678910111213141516171819202122232425262728293031323334353637383940414243444546474849505152535455565758596061626364656667686970717273747576777879808182838485868788899091929394959697989900

**K.** Download the User manual and consult the Intended purpose.

**N.** Units SI Imperial

**O.** The result of the algorithm's calculation is: ... points. Set all parameters to calculate prediction. Here a short section will be provided to help with the result interpretation. This piece of text can be general for all results, or can be shown depending when the certain conditions are met. This can include statement into which the risk classification the calculated result can be stratified (e.g. High, Moderate, Low). Also the performance data in the Internal and relevant External validation cohorts can be shown here such as but not limited to: the c-satistic, sensitivity, specificity together with the number of cases of the condition in scope within the cohort. Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer.

**Figure 1.** Example of an algorithm landing page on the Evidencio website.

**A. Algorithm title**

This is the title and name of the algorithm.

**B. Algorithm description**

This is a short description of the algorithm.

**C. Research authors**

These are the research authors of the paper that originally published the algorithm.

## D. Algorithm tags

These are the tags that are assigned to the algorithm. Evidencio has the following status tags: "Draft", "Public", "Private", "Under review". Evidencio has the following algorithm type tags: "Composite algorithm", "Sequential algorithm", "API algorithm". Evidencio has the following calculation method tags: "Linear regression", "Logistic regression", "Cox regression", "RScript" and "Custom calculation". Next to this, there are tags that indicate the specialty e.g. "Cardiology".

## E. LOT number

The LOT number indicated the algorithm version, the algorithm identifier, and the algorithm publication date. Publication date is indicated as YY.MM.DD.

Additionally, the CE mark is displayed next to the LOT number. This way, medical devices can be easily recognized.

## F. UDI Number

For information on the UDI Number see **Section 5.2** on **page 5** of this user manual.

## G. Details button

On the top right of the algorithm page, several clickable buttons are displayed that show a pop-up when clicked. The first button opens a pop-up concerning additional information about the algorithm. This pop-up has three sections: Details, Study characteristics and Supporting publications & related files.

### Details

The first part of the additional information concerns the details of the algorithm as shown in **Figure 2**. This section may show the calculation if it is built as a mathematical formula and, if applicable, shows the conditions at which certain formulas are used.

**Details**


Algorithm author	Evidencio	Status	Draft
Algorithm ID	10513	Share	
Version	1.0		
Revision date	2025-02-11		
Specialty	Cardiology , Geriatrics , Vascular medicine		
Algorithm type	Custom calculation (Conditional)		
MeSH terms	<ul style="list-style-type: none"> <li>• Heart Failure</li> <li>• Diabetes Mellitus</li> <li>• Elderly</li> </ul>		

Condition	Formula
Categorical Variable 1=Yes	$\text{Categorical Variable 1} + \text{Categorical Variable 2}^2 + \frac{3 \cdot \text{Continuous Variable 1}}{\text{Continuous Variable 2}}$
Categorical Variable 1=No	$\sqrt{\text{Continuous Variable 1}} + \frac{2 \cdot \text{Continuous Variable 2}}{\text{Continuous Variable 3}}$

**Figure 2.** Example of first part of Details section.

## Study Characteristics

Below the 'Details section' the section labelled "Study characteristics" provides information on the characteristics of the patient data used to derive and validate the algorithm. Additional information is provided on the methods used to develop and/or validate the algorithm. An example of the Study characteristics section can be seen in **Figure 3**.



**Figure 3.** Example of the study characteristics section under the Details tab.

## Supporting publications & Related files




An important part of the Study characteristics is the information on Supporting publications and related files. The list of related files and relevant tags can also be found in **Paragraph 8.5**. These sections can be found at the bottom of the Details-pop-up as shown in **Figure 4**.

### Supporting Publications

<p><b>Title or description</b> Title Derivation Paper DOI: DOI: 10.1234/ABCD.1234.5678</p> <p>Title External Validation DOI: DOI: 10.1234/ABCD.1234.5678</p> <p>Title Peer Review Paper DOI: DOI: 10.1234/ABCD.1234.5678</p>	<p><b>Tags</b></p> <ul style="list-style-type: none"> <li>Original calculator</li> <li>Internal validation</li> <li>External validation</li> <li>Peer review</li> </ul>
--	---

### Related files

<p><b>Preview</b></p>   	<table border="0"> <thead> <tr> <th style="text-align: left;">Name</th> <th style="text-align: left;">Tags</th> </tr> </thead> <tbody> <tr> <td>Derivation Paper.pdf 24.93 KB</td> <td>Original calculator Internal validation External validation</td> </tr> <tr> <td>External Validation.pdf 24.93 KB</td> <td>External validation</td> </tr> <tr> <td>Peer Review Paper.pdf 24.93 KB</td> <td>Peer review</td> </tr> </tbody> </table>	Name	Tags	Derivation Paper.pdf 24.93 KB	Original calculator Internal validation External validation	External Validation.pdf 24.93 KB	External validation	Peer Review Paper.pdf 24.93 KB	Peer review	<p><b>Tags</b></p> <ul style="list-style-type: none"> <li>Original calculator</li> <li>Internal validation</li> <li>External validation</li> <li>Peer review</li> </ul>
Name	Tags									
Derivation Paper.pdf 24.93 KB	Original calculator Internal validation External validation									
External Validation.pdf 24.93 KB	External validation									
Peer Review Paper.pdf 24.93 KB	Peer review									

**Figure 4.** Example of the Supporting publication & Related files section under the Details tab.

## H. Intended purpose

Under this tab, the intended purpose can be found, containing a lot of information regarding the algorithm, its user, target population, clinical benefit, etc. This information is also provided in this manual and can be found in **Chapter 6 on page 5**.

## I. Electronic label

The electronic label button opens a pop-up with the location and address of Evidencio, the LOT number, the UDI number, the CE-mark, the medical device logo and a download link for the declaration of conformity of the medical device. The example of the electronic label is shown in **Figure 5**. The electronic label is unique for each algorithm comprising the Child-Pugh Score.

### Extra Information

Intended purpose [Electronic label](#) [Release Notes](#)

---

Algorithm Name

 Evidencio B.V., Irenesingel 19, 7481 GJ Haaksbergen, The Netherlands

 V-1.0-10513.25.03.19

 (01)0123456787654321AB(8012)v1.0(4326)250319(240)10513

 Download the [User manual](#)

 In vitro diagnostic medical device

 Download the [Declaration of conformity](#)

**Figure 5.** Example of an electronic label under the Electronic Label tab.


## J. Release notes

Under this tab the most recent release notes can be found, noting the most significant changes between the versions of the algorithm found on the Evidencio website.


The 'Release Notes' button opens a pop-up with the latest release notes of the algorithm. Here you can find a list of the most significant changes over the different versions of the algorithm. Additionally, if there are any known residual anomalies the user should be aware of, they are listed here. It is recommended to read these notes after a version update to see if these changes are relevant to you.

## K. User manual

This user manual can be found in three places: 1) under the short description of the algorithm on the Evidencio algorithm page, 2) on the right of the algorithm page, and 3) as a tab in the electronic label screen. Additionally, all versions of the user manual can be found in the general page for all user manuals for medical devices. The page can be found under the 'About' drop-down menu button as shown in **Figure 6**. The user manual page is shown in **Figure 7**. This version of the manual can be printed if required. If necessary, a paper version of the manual can be requested to be sent to you by mail. Evidencio's contact details are listed in **Chapter 12** of this user manual.



**Figure 6.** The drop-down menu where the user manual page can be found.



**Figure 7.** The user manual page for all user manuals.

## L. Languages


Here an overview of languages in which the Child-Pugh Score is available is provided, any of which can be selected by clicking on the corresponding flag icon. The standard language on the Evidencio website is English.

Please note that, if a language is selected, only the user interface of the specific algorithm will be translated, other general features and information on the site might still be set to one of our primary languages English, German, and Dutch.

When you find mistranslations, irregularities, confusing or ambiguous use of language in English or any other language on the Evidencio website or in one of our manuals, please do not hesitate to contact us using the contact information provided at the end of this manual.

## M. Version selection

If available, clicking on the Version tab allows the user to select a different version of the Child-Pugh Score from a list as displayed in **Figure 8**. Please note that the algorithm currently selected is not presented in the dropdown menu.




**Figure 8.** Example of version selection tab.

## N. Input section


The Evidencio platform allows two separate input variables; categorical variables and continuous variables.

### Categorical variables

In the example shown in **Figure 9** and **Figure 10**, the example **Categorical Variable 1** concerns a categorical variable. The input that is wished to be used can be entered by clicking on either button. The selected button changes to green, as seen in **Figure 10**.



**Figure 9.** Example of a categorical variable, no button has been clicked and thus no input has been provided by the user.




**Figure 10.** Example of a categorical variable, where the "Yes" button has been clicked.

### Continuous variables

In the example shown in **Figure 11**, the **Continuous Variable 3**, exemplifies a continuous variable. The plausible ranges for which the algorithm is tested and deemed valid are used.


The details for a patient can be entered by sliding the button to the correct value, or by entering the correct value in the box on the right-hand side (i.e., where the 10.2 mg/dL is entered for the **Continuous Variable 3**).



**Figure 11.** Example of a continuous variable, where "10.2 mg/dL" has been entered.

## Unit conversion

Sometimes it is possible to use a unit conversion, by clicking on the unit when the green arrows are present. See **Figure 12** below where the unit has been clicked and switched.



**Figure 12.** Example of a continuous variable where “50.1  $\mu\text{mol/L}$ ” has been entered.

## Details on variable measurements

Directly underneath the name for each variable, additional details can be provided on, for example, the methods required to enter the correct value for each variable. Details may include but are not limited to; more detailed explanation of the variable, the ranges of the variables (for healthy individuals), or a description when a continuous variable should be true or false (cut-off values).

## O. Result section

At the bottom of the algorithm landing page, the results of the algorithm are shown.

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer on: <https://www.evidencio.com/disclaimer>.

## Result calculation

When all variables are filled in, and the user presses calculate, a result can be calculated. No result is displayed until all variables are filled in and the result section will indicate; “Set all parameters to calculate prediction.”.

## Result interpretation

In the result interpretation, a stratification may be provided based on the calculated results. Additional information about this stratification and the classification as found in the derivation and important validation cohorts may also be provided. An example of the information is shown in **Figure 13**.

**The result of the algorithm's calculation is: ... points.**

Set all parameters to calculate prediction.

Here a short section will be provided to help with the result interpretation. This piece of text can be general for all results, or can be shown depending when the certain conditions are met.

This can include statement into which the risk classification the calculated result can be stratified (e.g. **High, Moderate, Low**).

Also the performance data in the Internal and relevant External validation cohorts can be shown here such as but not limited to; the c-satistic, sensitivity, specificity together with the number of cases of the condition in scope within the cohort.

Calculations alone should never dictate patient care, and are no substitute for professional judgement. See our full disclaimer.

**Figure 13.** Example of the result display and information section.

## 10. Implementation of the algorithm through an API

The Child-Pugh Score can be used through Evidencio's API to allow for (automated) calculation of the severity of liver dysfunction. In the case of use of the MDSW through the API, the user should take into account the different inputs for the algorithm, in order to properly interpret the results.

The information provided over the API is the same as the information that is displayed in the graphical user interface on the web application provided by Evidencio. In **Box 1** below, an example of a result from the Child-Pugh Score over the API is shown. The result concerns a JSON formatted text. The API for the Child-Pugh Score leverages the generic API that is provided for the Evidencio platform and therefore contains information that may be applicable for different software algorithms and devices. This means that not all of the details provided over the API may be relevant for the Child-Pugh Score.

```
{
  "CIPercentage": 0,
  "id": 253,
  "author": "Evidencio",
  "title": "Child-Pugh Score",
  "variables": {
    "8120": 1,
    "7946": 40,
    "5959": 1,
    "9047336707": 1,
    "9856": 0,
    "4142277124": 6,
    "2647": 1
  },
  "min": 9,
  "max": 9,
  "additionalResultSet": [
    {
      "key": "<p>Child-Pugh class B</p>",
      "min": 7,
      "max": 9,
      "id": 4308799064
    }
  ],
  "mintxt": "9",
  "maxtxt": "9",
  "result": "9",
  "resultText": "The Child-Pugh score:",
  "postresultText": "Points",
  "formulaSegments": [],
  "conditionalResultArray": [],
  "conditionalResultText": "",
  "UDI": "(01)08720938015090(8012)v2.2(4326)260626(240)253",
  "medicalDevice": "This is an in vitro diagnostic medical device. The electronic label is available at:
https://www.evidencio.com/models/show/253?v=2.2,
  "userManual": "Always refer to the user manual for correct use of the in vitro diagnostic medical device.
  The user manual can be found at: https://www.evidencio.com/manuals"
}
```

**Box 1:** Example of an API output for the Child-Pugh Score.

**Table 6** shows a match between the separately listed items in the API output with the items listed on the graphical user interface on the Evidencio website (elaborated in **Chapter 9**).

**Table 6.** Agreement between the API output and the graphical user interface on the Evidencio website.

API Item	GUI item	Comment
CIPercentage	N/A	Not applicable for the Child-Pugh Score as this feature is not used for the Child-Pugh Score.
id	Algorithm ID under 'details' Id used in the URL (www.evidencio.com/models/show/10277)	The ID is the Evidencio specific identification number of the algorithm.
author	Algorithm author under 'details'	Name of the Evidencio user who created the algorithm on the Evidencio platform.
title	Title of the algorithm (part A of figure 1).	-
variables	Input variables and their entered value. (part N. of figure 1)	The API displays the variables as unique IDs.
min	N/A	Depicts the lowest value when the result of the algorithm is a range. Since the Child-Pugh Score always displays a single value as a result, this value is the same as the 'result'.
max	N/A	Depicts the highest value when the result of the algorithm is a range. Since the Child-Pugh Score always displays a single value as a result, this value is the same as the 'result'.
additionalResultSet	Includes the Child-Pugh Score classification within a specific interval.	e.g. "Result interval 7 to 9 Child-Pugh class B".
mintxt	N/A	Same as 'min', but as a string.
maxtxt	N/A	Same as 'max', but as a string.
result	The main result of the algorithm, the Child-Pugh Score.	-
resultText	The text displayed in front of the main result	e.g. "The Child-Pugh Score:"
postresultText	The text displayed behind the main result	e.g. "Points.
formulaSegments	N/A	N/A
conditionalResultArray	Result interpretation displayed beneath "Conditional Information" (section O. in figure 1).	The API result shows the raw HTML text that is rendered by the software used for the graphical user interface.
conditionalResultText	Result interpretation displayed beneath "Conditional Information" (section O. in figure 1).	This section is the same as ' <i>conditionalResultArray</i> ', but displayed as a single string.
UDI	Same as the UDI displayed in the GUI (section F. in figure 1).	-
medicalDevice	The electronic label (section I. in figure 1).	The API refers to the electronic label on the graphical user interface.
userManual	The user manual (section K. in figure 1).	The API refers to the location of the user manual at the user interface & Evidencio website).

Instructions on how to implement the API within a system are included in a separate document that is made available to the party performing the technical implementation. The party performing the integration of the Child-Pugh Score using the API should adhere to the requirements outlined in **253-DOC-45 Instructions for API integration Child-Pugh Score**.

NOTE: Specific input variables, such as Clotting time measurement, can be fixed for integrated use of Child-Pugh Score. In such cases, the user is unable to change the Clotting time measurement method (either "International Normalized Ratio (INR)" or "Prothrombin time").

## 11. User manual revision history

Version	Revision notes
V1.0	Original version
V2.0	Updated with newest versions after adding multiple languages Updated the manufacturer address

## 12. Manufacturer details

Contact details of Evidencio:



Evidencio B.V., Thermen 1, 7521 PS Enschede, The Netherlands

[www.evidencio.com](http://www.evidencio.com)

tel: +31 53 85195 08

e-mail: [info@evidencio.com](mailto:info@evidencio.com)